

ENDOACCESS

ENDO FOR LIFE



**ENDOSCOPIC ACCESSORIES**



## Disposable Instruments

<b>Disposable ESD Instruments</b>	<b>Page 1 - 3</b>
<b>Disposable ERCP Instruments</b>	<b>Page 3 - 5</b>
<b>Disposable EMR Instruments</b>	<b>Page 5 - 8</b>
<b>Disposable Cold Snares</b>	<b>Page 8</b>
<b>Disposable Snare with Bag</b>	<b>Page 9</b>
<b>Disposable Hemoclips</b>	<b>Page 9</b>
<b>Disposable Grasping Forceps</b>	<b>Page 10 - 13</b>
<b>Disposable Retrieval Baskets</b>	<b>Page 14 - 15</b>
<b>Disposable Biopsy Forceps</b>	<b>Page 15 - 19</b>
<b>Disposable Hot Biopsy Forceps</b>	<b>Page 19</b>
<b>Disposable Spray Catheters</b>	<b>Page 20</b>
<b>Disposable Cytology Brushes</b>	<b>Page 21</b>
<b>Disposable Cleaning Brushes</b>	<b>Page 21 - 22</b>
<b>Disposable Accessories</b>	<b>Page 23 - 24</b>

## Reusable Instruments

<b>Reusable Biopsy Forceps</b>	<b>Page 25 - 27</b>
<b>Reusable Foreign Body Removal</b>	<b>Page 28 - 31</b>
<b>Reusable Retrieval Baskets</b>	<b>Page 32 - 33</b>
<b>Reusable Electrical Snares</b>	<b>Page 33 - 35</b>
<b>Reusable Scissors Type Grasping Forceps</b>	<b>Page 36</b>
<b>Reusable Cleaning Brushes</b>	<b>Page 36</b>

## ESD Instruments

### Disposable Electrical Forceps (Rotatable, Hemo Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24HFR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24HFR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

mark, cut and stop bleeding by opening and closing the branches. Slow rotation facilitates operation  
Olympus compatible HF connection

### Disposable Electrical Forceps (Rotatable, Tooth/Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DTR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DTR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access. Hook tip designed for lifting of tissue  
Olympus compatible HF connection

### Disposable Electrical Forceps (Rotatable, Dissection)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DFR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DFR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access.  
Olympus compatible HF connection



## Disposable Electrical Forceps (Rotatable, Dissection Scissors)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DSR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DSR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access  
Olympus compatible HF connection

## Disposable Ball Type



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24BT18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	2mm
EA-S24BT23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	2mm

**Multifunctional Design:** marking, cutting, dissection, coagulation, spraying, etc

**Ball Shape Type:** Good operability, precise coagulation effect. Good traction during cutting and peeling

**Spray function:** For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

## Disposable Cylindrical Type (Tunnel Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24ST18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	1.5/2mm adjustable
EA-S24ST23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	1.5/2mm adjustable

**Multifunctional Design:** Marking, cutting, injection, lifting, dissection, coagulation, etc.

**Needle shape tip:** Sharp cutting tip extendable up to 2 mm

**Submucosa injection function:** High quality hydraulic solution pressure with needle point design; the knife can inject precisely into the submucosa and achieve a perfect mucosal lifting effect

Olympus compatible HF connection

## Disposable Ball Type (CT Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24B3T18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4mm
EA-S24B3T23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4mm

**Insulated Tip Design:** Ensures safety when cutting and dissecting, excellent conductivity; gentle cutting and stable coagulation

**Adjustable electrode length:** Adjustable from 0 mm to 4 mm according to operation procedure and lesions

**Spray function:** For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

## Disposable Cylindrical Knife (Square Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24S2T18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4mm
EA-S24S2T23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4mm

**Square Shape Tip Design:** Marking and coagulation

The design of the tip allows to lift or hook the tissue to cut efficiently. Good support surface and better coagulation effect


**Adjustable electrode length:** Adjustable cutting tip from 0 to 4mm

**Spray function:** For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

## Disposable ERCP Instruments

### Disposable Papillotome (rotatable)



Item	Tube Ø and Working Length	Working Channel Ø	Tip length A	Cutting length B
<b>Double Lumen</b>				
EA-S0720EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	20mm
EA-S0725EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	25mm
EA-S0730EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	30mm
<b>Triple Lumen</b>				
EA-S0720EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	20mm
EA-S0725EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	25mm
EA-S0730EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	30mm

braided wire, radiopaque markers ; preformed and tapered tip to Ø1.9mm

### Disposable Sphincterotome



Item	Tube Ø and Working Length	Working Channel Ø	Needle Length
EA-S24NE212 (Double Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	3-5mm (adjustable)
EA-S24NE213 (Triple Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	3-5mm (adjustable)

adjustable needle length 3-5mm, monofil, radiopaque marker

### Disposable ERCP Basket



Item	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24BS21-1 (Single Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	20mm
EA-S24BS21-2 (Double Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	20mm

For stone extraction, suitable for lithotripsy, hexagonal, retractable atraumatic tip, double lumen, basket length 45mm



## Disposable ERCP Cytology Brush



Item	Tube Ø and Working Length	Working Channel Ø
EA-S24CB21-2 double lumen, brush length 10mm	Ø2.4mm x 2100mm	≥Ø2.8mm

## Disposable Guide Wire



Item	Tube Ø and Working Length	Tip
EA-GW35450	Ø0.035" x 4500mm	55 mm hydrophil straight tip
EA-GW25450	Ø0.025" x 4500mm	55 mm hydrophil straight tip
EA-GW35450 J	Ø0.035" x 4500mm	55 mm hydrophil J tip
EA-GW25450 J	Ø0.025" x 4500mm	55 mm hydrophil J tip

## Disposable EMR Instruments

### Disposable Injection Needle




- 1) stainless steel needle point, teflon inner & outer tube
- 2) EA-SxxxxNJxx xx
  - | | | -----needle length
  - | | | ----- needle diameter
  - |----- tube diameter & working length

Item	Use	Tube Ø and Working Length	Working Channel Ø	Needle Length	Needle Ø
EA-S1810NJxxxx	Gastroscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	4/5mm	23/25G
EA-S1816NJxxxx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	4/5/6mm	23/25G
EA-S2416NJxxxx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2418NJxxxx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2421NJxxxx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2423NJxxxx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4/5/6mm	21/23/25G

with locking mechanism

## Disposable Electrical Snare (Oval Type) multifil rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
<b>Rotatable</b>				
EA-S18OS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18OS12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18OS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24OS16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm
<b>Non-rotatable</b>				
EA-S18OSN10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18OSN12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18OSN16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24OSN16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable

## Disposable Electrical Snare (Oval Type) monofil rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
<b>Rotatable</b>				
EA-S24OSM23xxR	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm
<b>Non-rotatable</b>				
EA-S24OSM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable




## Disposable Electrical Snare (Hexagonal Type) multifil rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
<b>Rotatable</b>				
EA-S18HS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HS12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24HS16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm
<b>Non-rotatable</b>				
EA-S18HSN10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HSN12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HSN16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24HSN16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection

## Disposable Electrical Snare (Half-Moon Type) multifil non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18HM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HM12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HM16xx	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	10/15/20mm
EA-S24HM16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/25/30/mm

Asymmetric, Olympus compatible HF connection

## Disposable Cold Snares

### Disposable Grasping Forceps (Cold Snare Type, monofil)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24CS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15mm

### Disposable Snare (Cold Snare Type with HF connection, monofil)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24OSC23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15mm

Olympus compatible HF connection

## Disposable Snare with Bag

### Disposable Snare with Bag (also available with net)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18SN10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	20mm
EA-S18SN12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	20mm
EA-S18SN16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	20mm
EA-S24SN18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	20/30mm
EA-S24SN23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/30mm

Tear proof, perforated extraction bag for removal of foreign bodies and polyps

## Disposable Hemoclip

### Disposable Hemoclip



Item	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24HC16xxR	Ø2.4mm x 1600mm	≥Ø2.8mm	9/11/14/17mm
EA-S24HC18xxR	Ø2.4mm x 1800mm	≥Ø2.8mm	9/11/14/17mm
EA-S24HC23xxR	Ø2.4mm x 2300mm	≥Ø2.8mm	9/11/14/17mm

Preloaded, repositionable, rotatable



## Disposable Grasping Forceps

### Disposable Grasping Forceps (Rat Tooth Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18RG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18RG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18RG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5mm
EA-S24RG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6mm
EA-S24RG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6mm
EA-S24RG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-S24RG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6mm

For the removal of flat and hard object (e.g. coins)

### Disposable Grasping Forceps (Rat Tooth with Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18RAG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18RAG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18RAG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5mm
EA-S24RAG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6mm
EA-S24RAG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6mm
EA-S24RAG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-S24RAG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6mm

Sharp teeth, long jaws for difficult retrieval of flat and hard objects

## Disposable Grasping Forceps (Pelican Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18PG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	9.5mm
EA-S18PG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	9.5mm
EA-S18PG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	9.5mm
EA-S24PG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	15mm
EA-S24PG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	15mm
EA-S24PG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	15mm
EA-S24PG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	15mm

For the retrieval of big objects or stent removal

## Disposable Grasping Forceps (Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10AG10	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	5.5mm
EA-S10AG12	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	5.5mm
EA-S18AG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	7mm
EA-S18AG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	7mm
EA-S18AG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	7mm
EA-S24AG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	7.5mm
EA-S24AG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7.5mm
EA-S24AG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	7.5mm
EA-S24AG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7.5mm

For the removal of sharp-edged objects

## Disposable Grasping Forceps (3-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18TP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18TP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18TP12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S24TP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24TP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24TP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24TP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies

## Disposable Grasping Forceps (3-Prong Hook Type)




Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18TP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18TP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18TP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

This item is also available on request in Ø 2.4. Longer lead times have to be considered.




## Disposable Grasping Forceps (5-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18FP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S24FP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24FP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24FP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24FP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies

## Disposable Grasping Forceps (5-Prong Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18FP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18FP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

This item is also available on request in Ø 2.4. Longer lead times have to be considered

## Disposable Grasping Forceps (2-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S24DP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24DP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24DP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24DP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

For the removal of gastric balloons

## Disposable Retrieval Baskets

### Disposable Retrieval Basket (Fold Angular, stainless steel)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18BS07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18BS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18BS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24BS15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

### Disposable Retrieval Basket (without Tip, Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10UB07xx	Urology	Ø1.0mm x 700mm	≥Ø1.2mm	10/15/mm
EA-S10UB10xx	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	10/15/mm
EA-S10UB12xx	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	10/15/mm
EA-S18BM07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18BM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18BM16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24BM15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

## Disposable Spirality Basket (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10U2B07xx	Urology	Ø1.0mm x 700mm	≥Ø1.2mm	10/15/mm
EA-S10U2B10xx	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	10/15/mm
EA-S10U2B12xx	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	10/15/mm
EA-S18B2M07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18B2M10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18B2M16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S24B2M15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

## Disposable Biopsy Forceps

### Disposable Biopsy Forceps (Oval, fenestrated cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10C10	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	4.5mm
EA-S10C12	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	4.5mm
EA-S18C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18C12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S18C16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5.5mm
EA-S24C16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24C18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm
EA-S30C23	Colonoscopy	Ø3.0mm x 2300mm	≥Ø3.2mm	6.5mm

Oval, fenestrated cups ensure minimal damage of tissue. Allows deeper biopsies. Uncoated



## Disposable Biopsy Forceps (Oval, coated, fenestrated cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S100C10	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	4.5mm
EA-S100C12	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	4.5mm
EA-S180C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S180C12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S180C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S180C16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5.5mm
EA-S180C18	Gastroscopy	Ø1.8mm x 1800mm	≥Ø2.0mm	5.5mm
EA-S240C16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S240C18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S240C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S240C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm
EA-S300C23	Colonoscopy	Ø3.0mm x 2300mm	≥Ø3.2mm	6.5mm

Oval, fenestrated cups ensure minimal damage of tissue. Allows deeper biopsies. Coated

## Disposable Biopsy Forceps (Oval Cups with Needle, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S180W10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S180W12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S180W15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S240W16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.0mm	6.5mm
EA-S240W18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S240W21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S240W23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping. Minimal damage of tissue. Uncoated

## Disposable Biopsy Forceps (Oval Cups with Needle, coated, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18CW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18CW12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18CW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24CW16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24CW18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24CW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24CW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping. Minimal damage of tissue. Coated

## Disposable Biopsy Forceps (Alligator Cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18A10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18A12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18A15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24A16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24A18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24A21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24A23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Uncoated

## Disposable Biopsy Forceps (Alligator Cups, coated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18AC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18AC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18AC15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24AC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24AC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24AC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24AC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Coated

## Disposable Biopsy Forceps (Alligator Cups with Needle, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18AW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18AW12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18AW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24AW16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24AW18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24AW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24AW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue. Uncoated



## Disposable Biopsy Forceps (Alligator cups with Needle, coated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18WA10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18WA12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18WA15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24WA16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24WA18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24WA21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24WA23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue. Coated

## Disposable Hot Biopsy Forceps

### Disposable Hot Biopsy Forceps



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18H10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18H12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18H15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-S24H16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24H18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24H21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24H23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

For the ablation of small lesions with simultaneous hemostasis or ablation of diminutive polyps. Olympus compatible HF connection

## Disposable Spray Catheters

### Disposable Standard Spray Catheter



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18SC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18SC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S18SC16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm
EA-S24SC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24SC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm
EA-S24SC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24SC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Even spray for clear images. For dye, saline and other solutions

### Disposable Fog Spray Catheters



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18FC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S18FC16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm
EA-S24FC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24FC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm
EA-S24FC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24FC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Extensive spray. For dye, saline and other solutions

## Disposable Cytology Brushes

### Disposable Cytology Brushes



Item	Use	Tube Ø and Working Length	Working Channel Ø	Ø Brushes and Length
EA-S18CB10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	3mm/7mm
EA-S18CB12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	3mm/7mm
EA-S18CB15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	3mm/7mm
EA-S24CB16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	3mm/7mm
EA-S24CB18	Colonoscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	3mm/7mm
EA-S24CB23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	3mm/7mm

Excellent brush hardness for improved sample collection. Thinner bristles for significantly improved cell collection capability.

## Disposable Cleaning Brushes

### Disposable Cleaning Brush (Double ended)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S10DB12	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm
EA-S10DB22	Urology	Ø1.0mm x 2200mm	≥Ø1.2mm
EA-S18DB18	Bronchoscopy	Ø1.8mm x 1800mm	≥Ø2.0mm
EA-S24DB22	Colonoscopy	Ø2.4mm x 2200mm	≥Ø2.8mm

For cleaning the working channel, symmetrical



## Disposable Cleaning Brush (Double ended, asymmetrical)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Ø Brushes
EA-S18DB23A35	Bronchoscopy	Ø1.8mm x 2300mm	≥Ø2.0mm	3mm/5mm
EA-S24DB22A510	Colonoscopy	Ø2.4mm x 2200mm	≥Ø2.8mm	5mm/10mm

For cleaning the working channel, asymmetrical

## Disposable Port Hole Cleaning Brush (Double ended, asymmetrical)

Valve cleaning



Item	Ø Brushes
EA-S1001PCD	5/10mm

## Disposable Port Hole Cleaning Brush

Valve cleaning



Item	Ø Brush
EA-S1001PC	10mm

## Disposable Accessories

### Disposable Colonoscopy Pants



Item	Size
EA-SCP-L	L (blue) (unit 200 pcs)
EA-SCP-XL	XL (green) (unit 120 pcs)

### Polyp Trap (4 chambers)



Item	Chambers
EA-SPT Removable chambers, transparent, rotatable	4
EA-SPT-S Transparent	1

### Disposable Helicobacter Test



Item
<b>EA-HEPY Urease Test</b>
Disposable helico check for the detection of H. pylori / results within 5 minutes

### Disposable Bite Block



Item	Strap
EA-SBBW	Yes
For adults	

### Disposable Bite Block (toothless adults)



<b>Item</b> EA-SBBK-N	<b>Strap</b> Yes
--------------------------	---------------------

### Disposable Bite Block (children)



<b>Item</b> EA-SBBK	<b>Strap</b> Yes
------------------------	---------------------

### Disposable Valve Cap



<b>Item</b> EA-SVC	<b>Compatibility</b> Olympus
-----------------------	---------------------------------

Closable lid



## Reusable Biopsy Forceps

### Reusable Biopsy Forceps (Standard Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18S07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5mm
EA-R18S10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-R18S15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-R24S15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24S17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24S21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24S23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

Standard, round fenestrated cups for minimal tissue damage, autoclavable

### Reusable Biopsy Forceps (Oval Fenestrated Cups Type)



Use	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18C07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24C15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24C17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups for deeper biopsy, autoclavable

## Reusable Biopsy Forceps (Oval Fenestrated Cups with Needle)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R180W07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5mm
EA-R180W10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-R180W15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-R240W15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	7.5mm
EA-R240W17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	7.5mm
EA-R240W21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	7.5mm
EA-R240W23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping.  
Minimal damage of tissue, autoclavable

## Reusable Biopsy Forceps (Alligator Jaws Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18A07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18A10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18A15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24A15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24A17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24A21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24A23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue, autoclavable

## Reusable Biopsy Forceps (Alligator, Single Moving Cup Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18OP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18OP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18OP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24OP15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24OP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24OP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24OP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Single tilting cup (90°) for tangential sampling, autoclavable

## Reusable Biopsy Forceps (Alligator Jaws with Needle)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18AW07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18AW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18AW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24AW15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24AW17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24AW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24AW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue. Uncoated, autoclavable



## Reusable Foreign Body Grasping Forceps

### Reusable Grasping Forceps (Pelican Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18PG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	11.5mm
EA-R18PG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	11.5mm
EA-R18PG15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	11.5mm
EA-R18PG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	11.5mm
EA-R24PG15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	15mm
EA-R24PG17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	15mm
EA-R24PG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	15mm
EA-R24PG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	15mm

For the retrieval of big objects or stent removal, autoclavable

### Reusable Grasping Forceps (Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18AG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	7mm
EA-R18AG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	7mm
EA-R18AG15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	7mm
EA-R18AG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	7mm
EA-R24AG15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	8mm
EA-R24AG17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	8mm
EA-R24AG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	8mm
EA-R24AG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	8mm

For the removal of sharp-edged objects, autoclavable

## Reusable Grasping Forceps (Rat Tooth Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18RG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RG15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24RG15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24RG17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24RG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24RG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

For the removal of flat and hard object (e.g. coins), autoclavable

## Reusable Grasping Forceps (Rat Tooth Type with Alligator Jaws Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RAG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18RAG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RAG15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24RAG15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24RAG17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24RAG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24RAG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

Sharp teeth, long jaws for difficult retrieval of flat and hard objects, autoclavable

## Reusable Grasping Forceps (3-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18TP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18TP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18TP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24TP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24TP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24TP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24TP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies, autoclavable


### Reusable Grasping Forceps (3-Prong Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18TP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18TP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18TP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

autoclavable

### Reusable Grasping Forceps (5-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18FP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18FP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18FP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24FP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24FP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24FP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24FP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies, autoclavable




## Reusable Grasping Forceps (5-Prong HookType)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18FP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18FP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18FP15H	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24FP23H autoclavable	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

## Reusable Grasping Forceps (2-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R24DP16	Gastrosocopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24DP17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24DP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24DP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

For the removal of gastric balloons, autoclavable

## Reusable Grasping Forceps (Rubber Tips Type)




Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RT10	Bronchoscopy	Ø 1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RT15	Gastrosocopy	Ø 1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24RT15	Gastrosocopy	Ø 2.4mm x 1500mm	≥Ø2.8mm	6mm
EA-R24RT17	Gastrosocopy	Ø 2.4mm x 1700mm	≥Ø2.8mm	6mm
EA-R24RT21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-R24RT23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm	6mm

For grasping slippery, flat and sharpedged objects, e.g. needles, wire, autoclavable

## Reusable Retrieval Baskets

### Reusable Fold Angular Basket (Stainless Steel)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18BS07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18BS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18BS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24BS15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable


### Reusable Basket without Tip (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18BM07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18BM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18BM16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24BM15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable

## Reusable Spirality Basket (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18SB07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18SB10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18SB16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24SB15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable

## Reusable Electrical Snares

### Reusable Electrical Snare (Oval Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18OS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18OS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24OS16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24OS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24OS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable, autoclavable

### Reusable Electrical Snare (Hexagonal Type)





Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18HS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18HS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24HS16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24HS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24HS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, autoclavable

## Reusable Electrical Snare (Half-Moon Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18HM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18HM16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R24HM16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	20/25/30mm
EA-R24HM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	20/25/30mm
EA-R24HM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/25/30mm

Asymmetric, Olympus compatible HF connection, autoclavable

## Reusable Scissors Type Grasping Forceps

### Reusable Grasping Forceps (Scissors Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R24ST15	Gastroscopy	Ø 2.4mm x 1500mm	≥Ø2.8mm	6mm
EA-R24ST17	Gastroscopy	Ø 2.4mm x 1700mm	≥Ø2.8mm	6mm
EA-R24ST21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-R24ST23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm	6mm

For cutting tissue, autoclavable

## Reusable Grasping Forceps (Half-Scissors Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R24HT15	Gastroscopy	Ø 2.4mm x 1500mm	≥Ø2.8mm
EA-R24HT17	Gastroscopy	Ø 2.4mm x 1700mm	≥Ø2.8mm
EA-R24HT21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm
EA-R24HT23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm

For cutting protruding tissue after ligature, autoclavable

## Reusable Cleaning Brushes

### Reusable Cleaning Brush



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18D16	Gastroscopy	Ø 1.8mm x 1600mm	≥Ø2.0mm
EA-R24D16	Colonoscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm
EA-R24D18	Colonoscopy	Ø 2.4mm x 1800mm	≥Ø2.8mm
EA-R24D22	Colonoscopy	Ø 2.4mm x 2200mm	≥Ø2.8mm

For cleaning the working channel,

### Reusable Porthole Cleaning Brush




Item
EA-R10PC


Valve cleaning

The logo consists of a yellow circle partially overlapping a black stylized shape that resembles a curved line or a partial 'G' shape.

**ENDOACCESS**  
ENDO FOR LIFE

A close-up photograph of a silver-colored metal surgical instrument with two serrated jaws, positioned in the upper left quadrant of the page.

**Endoaccess GmbH  
Feldriethe 1  
30826 Garbsen  
Germany**

A close-up photograph of a silver-colored metal surgical instrument, similar to the one in the top left, positioned in the bottom right corner of the page.

**+49 5131-44 22 61 0  
info@endoaccess.de  
www.endoaccess.de**